

# Care Gram

## Collection of Stool Sample

This following test(s) have been ordered by your doctor:

- C diff toxin – plain sterile container, **refrigerate** or deliver within 2 hours
- H pylori stool – plain sterile container, **refrigerate** or deliver within 2 hours
- Qualitative fecal fat – plain sterile container, **refrigerate** or deliver within 2 hours
- Culture (SSYCE) – Orange container – room temperature **\*\*Do not refrigerate\*\***
- O&P – Yellow SAF container, **room temperature or refrigerate**
- Giardia – Yellow SAF container, **room temperature or refrigerate**
- Cryptosporidium – Yellow SAF container, **room temperature or refrigerate**
- Other (fecal NA, K, Lytes) – plain sterile container, liquid specimen only, **refrigerate**
- Do not allow urine, toilet water or toilet tissue to mix with stool.
- Children – You may place plastic wrap inside the diaper and then transfer the stool to the container(s) provided. Do not bring the diaper with stool in it to the Laboratory.
- Do not collect stool if you have had a Barium enema in the last 14 days.

### Collection:

1. Place the collecting container (hat) under the seat and over the bowl of your toilet. The container should be dry.
2. Your bowel movement should go into the container.
3. Examine the stool for blood and mucus. Portions of stool with blood and mucus should be transferred to the specimen containers.
4. With the stick provided, add the stool in the plain container first.
5. The orange and yellow containers have spoons attached to the lids. Use the spoon to transfer the stool into the containers. Do not mix the lids. Add enough stool to bring the fluid up to the line.
6. Tighten the lid securely, and then mix the fluid and the stool by shaking the container.
7. Label the container(s) with your name, date of birth, collection date and time. Place the container(s) into the appropriate bag(s) provided and place the container(s) in the storage required (see above).
8. Wash your hands well with soap and water.
9. Deliver the specimen with the orders to the LMH Outpatient Laboratory area as soon as possible after the collection is over.

### Finding Out the Results

If you have any questions about your test, contact the LMH Laboratory at (740) 348-4155.

The results will be ready in 24 to 72 hours. You will need to check with the doctor who ordered the test(s) for the results.



**Licking Memorial  
Health Systems**

1320 West Main Street • Newark, OH 43055  
(740) 348-4000 • [www.LMHealth.org](http://www.LMHealth.org)

1616-0553  
09/20/2011